

Science Table

The Adjustable Height Science Table is ideal for incorporating computer technology in the classroom. This table has the ability to adapt easily to changing activities and technology by adding a wide variety of options.

Standard Features:

- Work surface height adjustment from 30" – 40" allows for multiple uses
- Develop your own arrangement to be compatible with the size and layout of your room
- Load bearing capacity of 500 lbs

Options:

- Height adjustment options - integral electric motor or may be powered with hand held cordless drill
- CPU holder suspends computer towers below work surface
- Hooks allow book bags to hang below work surface
- Modesty panels available in powder coated steel and wood veneer to fit customer requirements
- Dimensions may be adjusted to suit customer requirements



mott
MANUFACTURING

ISO 9001:2000 REGISTERED

QUALITY BY DESIGN

452 Hardy Rd, Brantford ON, Canada N3T 5L8 • P.O. Box 1120 N3T 5T3 • T. 519.752.7825 • F. 519.752.2895

www.mott.ca