

# Glassware Drying Unit

## Specialty Cabinet

These drying cabinets are used to dry and store glassware.

### Standard Features:

- Double wall insulated construction
- One adjustable stainless steel perforated shelf (not shown) for convenient loading / unloading. Additional shelves are available
- Stainless steel drip tray in cabinet bottom to contain spills (not shown)
- Pilot light
- 750 watt heating element with adjustable temperature up to 80°C / 176°F
- Removable control allows element replacement from front of cabinet
- Rear mounted element is located behind grill
- Electrical connections are located outside heated area
- This unit requires one 15 amp, 120 volt receptacle on the wall behind the cabinet
- Cabinet available in several different widths, depths and heights



**mott**  
MANUFACTURING

ISO 9001:2000 REGISTERED

QUALITY BY DESIGN

452 Hardy Rd, Brantford ON, Canada N3T 5L8 • P.O. Box 1120 N3T 5T3 • T. 519.752.7825 • F. 519.752.2895

[www.mott.ca](http://www.mott.ca)