

Mobile Sink Cabinet

Specialty Cabinet

This cabinet is designed for maximum versatility by providing a flexible supply of water in any research, school or care facility. This mobile sink solution can be used all through the facility when access to water is not available.

Hands Free Operation - Automatic faucet with infrared sensor provides hands-free operation.

Mobility - Plate mounted, non-marking 5" total lock swivel casters with santoprene rubber wheels allow the user to quickly and easily move to another work station, as tasks change.

Moisture & Corrosion Resistance - Available in either type 304 stainless steel or powder coated steel with protective Sigma Barrier™ coating. Durable stainless steel worktop with raised perimeter marine edge with integrally welded sink allows for easy cleaning.

Water Supply - 5 US gallon fresh water container. Also, 5 US gallon waste water container with drain spigot for easy release of waste water into the floor drain.

Operating Mode - The automatic pump provides up to 1.1 US gallons per minute with maximum current draw of 1.9 amperes at 12volt.

Battery & Charger - Onboard sealed 12volt battery and 1.5amp automatic charger including 6 inch cord and plug. Manual power switch on cabinet side.

Service Doors - Doors on front and back of cabinet provides convenient access for servicing.

Optional Equipment - Manual eyewash, manual faucet, supply hose to hook up with facility water supply.

* Available in various sizes, contact your local Mott dealer for more information.

** Stainless Steel cabinet shown.



Self Contained Water Supply



Water Delivery System



Infrared Sensor

mott
MANUFACTURING

ISO 9001:2000 REGISTERED

QUALITY BY DESIGN

452 Hardy Rd, Brantford ON, Canada N3T 5L8 • P.O. Box 1120 N3T 5T3 • T. 519.752.7825 • F. 519.752.2895

www.mott.ca