

Recessed Cabinets

Specialty Cabinet

Recessed Cabinets are typically found in hospital operating rooms. These storage cabinets maximize space utilization by safely storing the operating room inventory in an organized manner and serves as productive work areas.

Standard Features:

- Each cabinet can be constructed of type 304 or type 316 stainless steel with No. 4 finish. Cabinets are also available in powder coated steel.
- All exposed welds ground and polished to match No. 4 finish.
- Shelves are adjustable in 1/2" increments to deliver maximum adaptability.
- Available with full height hinged solid or framed glass doors with 6mm laminated safety glass. Lockable doors is an available option.
- Different sizes and configurations available, see our catalog for details.
- Stainless steel cabinets shown.



Custom desk/drawer cabinet features recessed task lighting and provisions for electrical outlets



mott
MANUFACTURING

ISO 9001:2000 REGISTERED

QUALITY BY DESIGN

452 Hardy Rd, Brantford ON, Canada N3T 5L8 • P.O. Box 1120 N3T 5T3 • T. 519.752.7825 • F. 519.752.2895

www.mott.ca